

MolecularDevices Nano

- Objectives
- Filter specifications
- Plate calibration
- Starting the system
 - Troubleshooting network problems
- Shutting down the system

Objectives

- 4x
- 10x (pixel size 0.655 μm)
- 20x

Filter specifications

Filter	Wavelengths	Shading correction plate
DAPI	Excitation: 377/50 nm	Green
	Emission: 447/60 nm	Pink
	Dichroic: 409 nm	
GFP	Excitation: 472/30 nm	Pink
	Emission: 520/35 nm	Red
	Dichroic: 495 nm	
Texas Red	Excitation: 562/40 nm	Red
	Emission: 624/40 nm	
	Dichroic: 593 nm	
Cy5	Excitation: 628/40 nm	Red
	Emission: 692/40 nm	
	Dichroic: 660 nm	

Plate calibration

Plate calibration is done in two steps:

1. Measure variation of plate bottom thickness
 - requires an empty plate
2. Measure required exposure time for autofocus laser
 - autofocus is based on the reflection of the laser from the plate and well bottom
 - reflection from the well bottom depends on the refractive index of the mounting media
 - requires a plate with relevant sample

The system has now these calibrated settings:

Manufacturer	Plate model	Objective and mounting
4titude	VisionPlate 96-well	
Costar	3596 96-well	10x PBS
Falcon	353072 96-well	10x, water
Greiner	Cellstar 655 180 96-well	20x PBS

Nunc	Optical 96-well	10x, air
PerkinElmer	CellCarrier Ultra 96-well	
PerkinElmer	ViewPlate 96-well	10x, PBS
PerkinElmer	CellCarrier Ultra 384-well	10x, PBS

If you use any other plates, you need to bring two plates for the training:

1. empty plate
2. a plate with sample (and media in the wells with no sample)

Starting the system

1. Power on the instrument.
2. Start computer.
3. Log in with account "Nano"
4. Mount L: disk
 - a. Double-click "Imu-L-disk" on the Desktop.
 - b. In the command window, write your university username in the format ATTK\username, press Enter.
 - c. In the command window, write your university password, press Enter.
5. Start MetaXpress
 - a. Double-click "MX" icon on the Desktop.
 - b. Select your data source (named after your university username).

Troubleshooting network problems

1. Check with LMU staff that you have correctly mounted the L-disk.
 - a. If L-disk is not available, use the local F disk (selector is in the "Run" tab).
2. Check if you can connect to your usual data source
 - a. If your usual data source doesn't work, select data source **zzz_<username>_local**, which uses the local database.

Shutting down the system

1. Close software
2. Disconnect L: disk (to delete login credentials)
 - a. Start Menu Computer
 - b. Right-click Network Location Imu (L:)
 - c. Select Disconnect
3. Shut down computer
4. Power off the instrument