

Blue Nuclear Label Red-shift

Sometimes the bleed-through of the blue nuclear label into the green channel seems higher than usual. There is a real phenomenon behind this increase; it's called photoconversion.

It is known for the blue nuclear labels DAPI, Hoechst and Vybrant DyeCyle Violet, but can happen to other labels too. Red-shifts of up to 80 nm have been observed in real samples, e.g. cy5-Beta-catenin on polyvinylidene (PVDF) membrane.

Note that the in-solution values stated in reference may differ from spectra in other conditions.

DAPI, Hoechst 33258 and Vybrant DyeCyle Violet Exhibit Photoconversion

- Photoconversion by exposure to UV light
 - The 405 nm laser is far enough from the UV to not cause this
- Photoconverted derivatives have low photostability
 - Bleached easily
 - Photobleaching of DAPI & Hoechst 33258 recovers 50% in 60 min
 - Photoconverted derivatives bleach at the same rate
- The photoconversion is reversible
 - After the photoconversion reversal an equilibrium state lasts for at least 2 more hours
- The photoconversion is dose dependent
 - The dose response of DAPI is stronger
 - 4x increase in 100 s
 - (Hoechst 2x in the same time)
- Hydroxyl peroxide induces similar photo-oxidation
 - Red-shift
 - Takes 4 h for Hoechst
 - 1.5 h for DAPI
 - Reversible

Spectral Characteristics of the Red-shifted Label

- Excitation at 450-490 nm
- Emission at 480-600 nm
 - Maxima at 540 nm

Too High Label Concentrations Can Diminish Signal

- High concentrations of DAPI & Hoechst 33258 diminish blue signal
- Self-quenching
- Signal in RNA increased

DAPI Spectral Shift

- DAPI binding with DNA is solution specific
- It can change from groove-binding to intercalation
- It can also induce polymer-dye adducts
- Emission maximum blue-shifts when complexed with DNA

DAPI Spectral Characteristics

- Excitation maxima at 364 nm
- Emission maximum at 454 nm

DAPI Red-shift Conditions

- DAPI also binds to RNA
- Emission maximum red-shifts when bound to RNA
 - Binding to RNA minor groove
 - Complexation with polyA-polyU
- Dissociation 100 times faster than from DNA
- High dye/phosphate ratio (1:10) reveals a new binding form
- New emission maximum 540 nm
- Binding of DAPI to sites in polynucleotides in proximity to previously bound DAPI
 - Disappears if background electrolyte is 0.4 M KCl
 - Disrupts DAPI-DAPI electrostatic interactions
 - DAPI-polyadenylic acid (polyA)
 - DAPI-Polyphosphate (polyP)
 - Excitation 415 nm
 - Emission 550 nm
 - Inositol phosphates (IP5 & IP6)
 - Emission 550 nm
 - Heparin

- Emission 550 nm
- Amorphous calcium phosphate (ACP)
 - Emission 550 nm
 - Not Hydroxyapatite (HAp)
- "Red-shifted DAPI fluorescence is not due to specific substrate chemistry, but indicates the presence of a high density of negatively-charged surfaces or molecules that locally concentrate DAPI. This local, increased DAPI concentration enables DAPI–DAPI interactions and its resultant red-shifted fluorescence."
- High concentrations of glycerol (as mounting media) are proportional to the level of DAPI photoconversion
- DAPI photoconversion rate has two components
 - Fast half-life < 10 s
 - Slow half-life > 60 s
 - Both are UV illumination intensity dependent (non-linear)

Hoechst 33258 Red-shift

- Exposure to UV causes both photobleaching and photoconversion
 - New emission maxima 540 nm
- Protonation of the dye
 - Also by exposure to hydrogen peroxide
- Similar properties in acidic environments (pH 0.5-3.0)
 - Quantum yield drops 80-fold
- QY increased 20-fold in pH 4.5
- Acid treatment reversible
- Equilibrium at 60 min past UV exposure
- Red-shifted protonated form shows up in nucleoli
 - RNA-binding
 - Excitation 369 nm
 - Emission 437 nm

Hoechst Spectral Characteristics

- Excitation 355 nm
- Emission 465 nm

Vybrant DyeCyle Violet

Vybrant DyeCyle Violet also exhibits the photoconversion red-shift of its spectra.

Mitigating the Red-shift

- Minimize UV intensity when imaging the blue nuclear labels.
 - Total dose or exposure time
- Use low label concentration
- Use low glycerol concentration in mounting media.
- Acquire blue nuclear channels last to minimize green channel false positives
- Use alternate nuclear dyes like DRAQ5 or RedDot which fluoresce in the far red region

Sources

- "A cautionary (spectral) tail: red-shifted fluorescence by DAPI–DAPI interactions"; Sidney Omelon, John Georgiou, Wouter Habraken; *Biochemical Society Transactions* Feb 09, 2016, 44 (1) 46-49; DOI: 10.1042/BST20150231
- "UV-induced Spectral Shift and Protonation of DNA Fluorescent Dye Hoechst 33258"; Dominika urek-Biesiada & Piotr Waligórski & Jurek W. Dobrucki; *J Fluoresc* (2014) 24:1791–1801; DOI 10.1007/s10895-014-1468-y
- "UV-Activated Conversion of Hoechst 33258, DAPI, and Vybrant DyeCycle Fluorescent Dyes into Blue-Excited, Green-Emitting Protonated Forms"; Dominika Zurek-Biesiada, Sylwia Kedracka-Krok, Jurek W. Dobrucki; *Cytometry Part A* (2013) 83A: 441-451; DOI: 10.1002/cyto.a.22260
- "High UV Excitation Intensity Induces Photoconversion of DAPI During Wide-Field Microscopy"; Stefan Rodic, Claire Brown, Erika (Tse-Luen) Wee; *McGill Science Undergraduate Research Journal* 9.1 (2014)
- Tanious, Farial A., et al. "DAPI (4', 6-diamidino-2-phenylindole) binds differently to DNA and RNA: minor-groove binding at AT sites and intercalation at AU sites." *Biochemistry* 31.12 (1992): 3103-3112.
- "Inositol phosphates induce DAPI fluorescence shift."; Kolozsvari, Parisi, Saiardi; *Biochem J.* 2014 Jun 15;460(3):377-85.; doi: 10.1042/BJ20140237.
- "Chromatic shifts in the fluorescence emitted by murine thymocytes stained with Hoechst 33342"; Timothy W. Petersen, Sherrif F. Ibrahim, Alan H. Diercks, and Ger van den Engh; *Cytometry Part A* (2004): 173-181; DOI: 10.1002/cyto.a.20058